



optiLube active CHGFree™

Meeting patient needs.
Responding to clinician demand.
Giving you a choice.

Optilube Active CHG Free contains local anaesthetic (Lidocaine 2%). It's ideally suited to patients with allergies to Chlorhexidine Gluconate or where its use is contraindicated.



6ml 11ml

Provided in a familiar syringe format
for urethral and suprapubic application

- Syringe falls easily from its protective tray, supporting an Aseptic Non Touch Technique when preparing a sterile field
- Manufacturing, expiry dates, and LOT numbers, on each syringe's tear off strip are easy to retain for accurate recording
- Easy-to-prime syringe for precise, controlled application
- Ergonomic design of the tip eliminates the risk of suction or damage to tissue
- Measuring facility on side of the syringe enables precise application of a controlled volume of jelly
- Latex-free and hypoallergenic

OM Code	Description	Inner pack	Outer pack
1167	Optilube Active CHG Free 6ml pre-filled syringe	Pack of 10	Case of 200
1168	Optilube Active CHG Free 11ml pre-filled syringe	Pack of 10	Case of 200

CE 0473
STERILE R



optimummedical®

Optimum Medical

First floor, Unit 5, Armley Court, Armley Road, Leeds, LS12 2LB, UK
Tel.: +44 (0) 8456 435 479, Fax: +44 (0) 1142 383 826
enquiries@optimummedical.co.uk • www.optimummedical.co.uk